

## Antibodies Product Datasheet

|                               |                                                                                                                                     |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | LIMK-1/2 (Phospho-Thr508/505) Antibody                                                                                              |
| <b>Catalog Number:</b>        | A0047                                                                                                                               |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                   |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                              |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/10000. Not yet tested in other applications.         |
| <b>Modification:</b>          | Phospho Specific                                                                                                                    |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.               |
| <b>Immunogen:</b>             | Synthesized phospho-peptide around the phosphorylation site of human LIMK-1/2 (phospho Thr508/505)                                  |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                             |
| <b>Gene Name:</b>             | LIMK1/LIMK2                                                                                                                         |
| <b>Protein Name:</b>          | LIM domain kinase 1/LIM domain kinase 2                                                                                             |
| <b>Other Name:</b>            | LIMK1; LIMK; LIM domain kinase 1; LIMK-1; LIMK2; LIM domain kinase 2; LIMK-2                                                        |
| <b>Human Swiss Prot No.:</b>  | P53667/P53671                                                                                                                       |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                             |
| <b>Specificity:</b>           | Phospho-LIMK-1/2 (T508/505) Polyclonal Antibody detects endogenous levels of LIMK-1/2 protein only when phosphorylated at T508/505. |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                        |