

Antibodies Product Datasheet

Product Name:	KOR-1 (Phospho-Ser369) Antibody
Catalog Number:	A0448
Reactivity:	Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from mouse KOR-1 around the phosphorylation site of Ser369. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	OPRK1
Protein Name:	Kappa-type opioid receptor
Other Name:	OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-KOR-1 (S369) Polyclonal Antibody detects endogenous levels of KOR-1 protein only when phosphorylated at S369.
Storage and Stability:	-20°C/1 year