

Antibodies Product Datasheet

Product Name:	CDK4 (Phospho-Thr172) Antibody
Catalog Number:	A0470
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CDK4 (Phospho-Thr172)
Concentration:	1 mg/ml
Gene Name:	CDK4
Protein Name:	CDK4 (Phospho-Thr172)
Other Name:	Cyclin-dependent kinase 4 (EC 2.7.11.22) (Cell division protein kinase 4) (PSK-J3)
Human Swiss Prot No.:	P11802
Formulation:	Liquid in PBS containing 50% glycerol, and 0.90% sodium azide.
Specificity:	This antibody detects endogenous levels of CDK4 (Phospho-Thr172) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year