

Antibodies Product Datasheet

Product Name:	ERF (Phospho-Thr526) Antibody
Catalog Number:	A0477
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ERF around the phosphorylation site of Thr526. AA range:492-541
Concentration:	1 mg/ml
Gene Name:	ERF
Protein Name:	ETS domain-containing transcription factor ERF
Other Name:	ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2
Human Swiss Prot No.:	P50548
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-ERF (T526) Polyclonal Antibody detects endogenous levels of ERF protein only when phosphorylated at T526.
Storage and Stability:	-20°C/1 year