

Antibodies Product Datasheet

Product Name:	CD130/gp130 (Phospho-Ser782) Antibody
Catalog Number:	A0838
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD130/gp130 (Phospho-Ser782)
Concentration:	1 mg/ml
Gene Name:	IL6ST
Protein Name:	CD130/gp130 (Phospho-Ser782)
Other Name:	Interleukin-6 receptor subunit beta (IL-6 receptor subunit beta) (IL-6R subunit beta) (IL-6R-beta) (IL-6RB) (CDw130) (Interleukin-6 signal transducer) (Membrane glycoprotein 130) (gp130) (Oncostatin-M receptor subunit alpha) (CD antigen CD130)
Human Swiss Prot No.:	P40189
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of CD130/gp130 (Phospho-Ser782) at Human:S782, Mouse:S780
Storage and Stability:	-20°C/1 year