

Antibodies Product Datasheet

Product Name:	CD150 (Phospho-Tyr327) Antibody
Catalog Number:	A0840
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD150 (Phospho-Tyr327)
Concentration:	1 mg/ml
Gene Name:	SLAMF1 SLAM
Protein Name:	CD150 (Phospho-Tyr327)
Other Name:	Signaling lymphocytic activation molecule (CDw150) (IPO-3) (CD antigen CD150)
Human Swiss Prot No.:	Q13291
Formulation:	Liquid in PBS containing 50% glycerol, and 0.191% sodium azide.
Specificity:	This antibody detects endogenous levels of CD150 (Phospho-Tyr327) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year