

Antibodies Product Datasheet

Product Name:	GATA-3 (Phospho-Ser308) Antibody
Catalog Number:	A0933
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GATA3 around the phosphorylation site of Ser308. AA range:274-323
Concentration:	1 mg/ml
Gene Name:	GATA3
Protein Name:	Trans-acting T-cell-specific transcription factor GATA-3
Other Name:	GATA3; Trans-acting T-cell-specific transcription factor GATA-3; GATA-binding factor 3
Human Swiss Prot No.:	P23771
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-GATA-3 (S308) Polyclonal Antibody detects endogenous levels of GATA-3 protein only when phosphorylated at S308.
Storage and Stability:	-20°C/1 year