

Antibodies Product Datasheet

Product Name:	HER2 (Phospho-Tyr1139) Antibody
Catalog Number:	A0943
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HER2 (Phospho-Tyr1139)
Concentration:	1 mg/ml
Gene Name:	ERBB2 HER2 MLN19 NEU NGL
Protein Name:	HER2 (Phospho-Tyr1139)
Other Name:	Receptor tyrosine-protein kinase erbB-2 (EC 2.7.10.1) (Metastatic lymph node gene 19 protein) (MLN 19) (Proto-oncogene Neu) (Proto-oncogene c-ErbB-2) (Tyrosine kinase-type cell surface receptor HER2) (p185erbB2) (CD antigen CD340)
Human Swiss Prot No.:	P04626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of HER2 (Phospho-Tyr1139) at Human:Y1139, Mouse:Y1140, Rat:Y1141
Storage and Stability:	-20°C/1 year