

Antibodies Product Datasheet

Product Name:	MAP-4 (Phospho-Ser696) Antibody
Catalog Number:	A1092
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAP4 around the phosphorylation site of Ser696. AA range:662-711
Concentration:	1 mg/ml
Gene Name:	MAP4
Protein Name:	Microtubule-associated protein 4
Other Name:	MAP4; Microtubule-associated protein 4; MAP-4
Human Swiss Prot No.:	P27816
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-MAP-4 (S696) Polyclonal Antibody detects endogenous levels of MAP-4 protein only when phosphorylated at S696.
Storage and Stability:	-20°C/1 year