

Antibodies Product Datasheet

Product Name:	L-plastin (Phospho-Ser5) Antibody
Catalog Number:	A1152
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human L-plastin (Phospho-Ser5)
Concentration:	1 mg/ml
Gene Name:	LCP1 PLS2
Protein Name:	L-plastin (Phospho-Ser5)
Other Name:	Plastin-2 (L-plastin) (LC64P) (Lymphocyte cytosolic protein 1) (LCP-1)
Human Swiss Prot No.:	P13796
Formulation:	Liquid in PBS containing 50% glycerol, and 0.108% sodium azide.
Specificity:	This antibody detects endogenous levels of L-plastin (Phospho-Ser5) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year