

## Antibodies Product Datasheet

|                               |                                                                                                                                          |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Cofilin (Phospho-Ser3) Antibody                                                                                                          |
| <b>Catalog Number:</b>        | A7047                                                                                                                                    |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                        |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                                   |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/10000. Not yet tested in other applications.              |
| <b>Modification:</b>          | Phospho Specific                                                                                                                         |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                        |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                    |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human Cofilin around the phosphorylation site of Ser3. AA range:1-50 |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                  |
| <b>Gene Name:</b>             | CFL1                                                                                                                                     |
| <b>Protein Name:</b>          | Cofilin-1                                                                                                                                |
| <b>Other Name:</b>            | CFL1; CFL; Cofilin-1; 18 kDa phosphoprotein; p18; Cofilin; non-muscle isoform                                                            |
| <b>Human Swiss Prot No.:</b>  | P23528                                                                                                                                   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                  |
| <b>Specificity:</b>           | Phospho-Cofilin (S3) Polyclonal Antibody detects endogenous levels of Cofilin protein only when phosphorylated at S3.                    |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                             |