

Antibodies Product Datasheet

Product Name:	Ribosomal Protein S6 (Phospho-Ser235) Antibody
Catalog Number:	A7214
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human S6 Ribosomal Protein around the phosphorylation site of Ser235. AA range:200-249
Concentration:	1 mg/ml
Gene Name:	RPS6
Protein Name:	40S ribosomal protein S6
Other Name:	RPS6; OK/SW-cl.2; 40S ribosomal protein S6; Phosphoprotein NP33
Human Swiss Prot No.:	P62753
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Ribosomal Protein S6 (S235) Polyclonal Antibody detects endogenous levels of Ribosomal Protein S6 protein only when phosphorylated at S235.
Storage and Stability:	-20°C/1 year