

Antibodies Product Datasheet

Product Name:	BLK (Phospho-Tyr389) Antibody
Catalog Number:	A8080
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BLK (Phospho-Tyr389)
Concentration:	1 mg/ml
Gene Name:	BLK
Protein Name:	BLK (Phospho-Tyr389)
Other Name:	Tyrosine-protein kinase Blk (EC 2.7.10.2) (B lymphocyte kinase) (p55-Blk)
Human Swiss Prot No.:	P51451
Formulation:	Liquid in PBS containing 50% glycerol, and 0.122% sodium azide.
Specificity:	This antibody detects endogenous levels of BLK (Phospho-Tyr389) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year