

Antibodies Product Datasheet

Product Name:	CAMKK1/2 (Phospho-Ser458/495) Antibody
Catalog Number:	A8086
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAMKK1/2 (Phospho-Ser458/495)
Concentration:	1 mg/ml
Gene Name:	CAMKK1 CAMKKA
Protein Name:	CAMKK1/2 (Phospho-Ser458/495)
Other Name:	Calcium/calmodulin-dependent protein kinase kinase 1 (CaM-KK 1) (CaM-kinase kinase 1) (CaMKK 1) (EC 2.7.11.17) (CaM-kinase IV kinase) (Calcium/calmodulin-dependent protein kinase kinase alpha) (CaM-KK alpha) (CaM-kinase kinase alpha) (CaMKK alpha)
Human Swiss Prot No.:	Q8N5S9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.150% sodium azide.
Specificity:	This antibody detects endogenous levels of CAMKK1/2 (Phospho-Ser458/495) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year