

Antibodies Product Datasheet

Product Name:	HIPK1/2/3 (Phospho-Tyr352/361/359) Antibody
Catalog Number:	A8111
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HIPK1/2/3 (Phospho-Tyr352/361/359)
Concentration:	1 mg/ml
Gene Name:	HIPK1 KIAA0630 MYAK NBAK2
Protein Name:	HIPK1/2/3 (Phospho-Tyr352/361/359)
Other Name:	Homeodomain-interacting protein kinase 1 (EC 2.7.11.1) (Nuclear body-associated kinase 2)
Human Swiss Prot No.:	Q86Z02
Formulation:	Liquid in PBS containing 50% glycerol, and 0.172% sodium azide.
Specificity:	This antibody detects endogenous levels of HIPK1/2/3 (Phospho-Tyr352/361/359) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year