

Antibodies Product Datasheet

Product Name:	MOK (Phospho-Tyr161) Antibody
Catalog Number:	A8146
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MOK (Phospho-Tyr161)
Concentration:	1 mg/ml
Gene Name:	MOK RAGE RAGE1
Protein Name:	MOK (Phospho-Tyr161)
Other Name:	MAPK/MAK/MRK overlapping kinase (EC 2.7.11.22) (MOK protein kinase) (Renal tumor antigen 1) (RAGE-1)
Human Swiss Prot No.:	Q9UQ07
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of MOK (Phospho-Tyr161) at Human:Y161, Mouse:Y161, Rat:Y161
Storage and Stability:	-20°C/1 year