

Antibodies Product Datasheet

Product Name:	NDR1/2 (Phospho-Thr444/442) Antibody
Catalog Number:	A8152
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDR1/2 (Phospho-Thr444/442)
Concentration:	1 mg/ml
Gene Name:	STK38 NDR1
Protein Name:	NDR1/2 (Phospho-Thr444/442)
Other Name:	Serine/threonine-protein kinase 38 (EC 2.7.11.1) (NDR1 protein kinase) (Nuclear Dbf2-related kinase 1)
Human Swiss Prot No.:	Q15208
Formulation:	Liquid in PBS containing 50% glycerol, and 0.176% sodium azide.
Specificity:	This antibody detects endogenous levels of NDR1/2 (Phospho-Thr444/442) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year