

Antibodies Product Datasheet

Product Name:	eIF2B-ε (Phospho-Ser540) Antibody
Catalog Number:	A8214
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human eIF2B-ε (Phospho-Ser540)
Concentration:	1 mg/ml
Gene Name:	EIF2B5 EIF2BE
Protein Name:	eIF2B-ε (Phospho-Ser540)
Other Name:	Translation initiation factor eIF-2B subunit epsilon (eIF-2B GDP-GTP exchange factor subunit epsilon)
Human Swiss Prot No.:	Q13144
Formulation:	Liquid in PBS containing 50% glycerol, and 0.159% sodium azide.
Specificity:	This antibody detects endogenous levels of eIF2B-ε (Phospho-Ser540) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year