

Antibodies Product Datasheet

Product Name:	Flt-1 (Phospho-Tyr1213) Antibody
Catalog Number:	A8228
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Flt-1 (phospho Tyr1213)
Concentration:	1 mg/ml
Gene Name:	FLT1
Protein Name:	Vascular endothelial growth factor receptor 1
Other Name:	FLT1; FLT; FRT; VEGFR1; Vascular endothelial growth factor receptor 1; VEGFR-1; Fms-like tyrosine kinase 1; FLT-1; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; FLT; Vascular permeability factor receptor
Human Swiss Prot No.:	P17948
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Flt-1 (Y1213) Polyclonal Antibody detects endogenous levels of Flt-1 protein only when phosphorylated at Y1213.
Storage and Stability:	-20°C/1 year