

Antibodies Product Datasheet

Product Name:	Gab 1 (Phospho-Tyr307) Antibody
Catalog Number:	A8245
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Gab 1 (phospho Tyr307)
Concentration:	1 mg/ml
Gene Name:	GAB1
Protein Name:	GRB2-associated-binding protein 1
Other Name:	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Human Swiss Prot No.:	Q13480
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Gab 1 (Y307) Polyclonal Antibody detects endogenous levels of Gab 1 protein only when phosphorylated at Y307.
Storage and Stability:	-20°C/1 year