

Antibodies Product Datasheet

Product Name:	Caveolin-2 (Phospho-Tyr27) Antibody
Catalog Number:	A8264
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caveolin 2 around the phosphorylation site of Tyr27. AA range:12-61
Concentration:	1 mg/ml
Gene Name:	CAV2
Protein Name:	Caveolin-2
Other Name:	CAV2; Caveolin-2
Human Swiss Prot No.:	P51636
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Caveolin-2 (Y27) Polyclonal Antibody detects endogenous levels of Caveolin-2 protein only when phosphorylated at Y27.
Storage and Stability:	-20°C/1 year