

Antibodies Product Datasheet

Product Name:	FANCG (Phospho-Ser383) Antibody
Catalog Number:	A8269
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human FANCG (phospho Ser383)
Concentration:	1 mg/ml
Gene Name:	FANCG
Protein Name:	Fanconi anemia group G protein
Other Name:	FANCG; XRCC9; Fanconi anemia group G protein; Protein FACG; DNA repair protein XRCC9
Human Swiss Prot No.:	O15287
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-FANCG (S383) Polyclonal Antibody detects endogenous levels of FANCG protein only when phosphorylated at S383.
Storage and Stability:	-20°C/1 year