

Antibodies Product Datasheet

Product Name:	GRK2 (Phospho-Ser685) Antibody
Catalog Number:	A8278sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho derived from human GRK2 (Phospho-Ser685)
Concentration:	1 mg/ml
Gene Name:	ADRBK1 BARK BARK1 GRK2
Protein Name:	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein coupled receptor kinase 2)
Other Name:	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein coupled receptor kinase 2)
Human Swiss Prot No.:	P25098
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This detects endogenous levels of GRK2 (Phospho-Ser685)
Storage and Stability:	-20°C/1 year