

Antibodies Product Datasheet

Product Name:	Ataxin-3 (Phospho-Ser256) Antibody
Catalog Number:	A8293
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Ataxin-3 (Phospho-Ser256)
Concentration:	1 mg/ml
Gene Name:	ATXN3 ATX3 MJD MJD1 SCA3
Protein Name:	Ataxin-3 (Phospho-Ser256)
Other Name:	Ataxin-3 (EC 3.4.19.12) (Machado-Joseph disease protein 1) (Spinocerebellar ataxia type 3 protein)
Human Swiss Prot No.:	P54252
Formulation:	Liquid in PBS containing 50% glycerol, and 0.193% sodium azide.
Specificity:	This antibody detects endogenous levels of Ataxin-3 (Phospho-Ser256) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year