

Antibodies Product Datasheet

Product Name:	BAG3 (Phospho-Tyr457) Antibody
Catalog Number:	A8296
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BAG3 (Phospho-Tyr457)
Concentration:	1 mg/ml
Gene Name:	BAG3 BIS
Protein Name:	BAG3 (Phospho-Tyr457)
Other Name:	BAG family molecular chaperone regulator 3 (BAG-3) (Bcl-2-associated athanogene 3) (Bcl-2-binding protein Bis) (Docking protein CAIR-1)
Human Swiss Prot No.:	O95817
Formulation:	Liquid in PBS containing 50% glycerol, and 0.170% sodium azide.
Specificity:	This antibody detects endogenous levels of BAG3 (Phospho-Tyr457) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year