

Antibodies Product Datasheet

Product Name:	BCLAF1 (Phospho-Ser531) Antibody
Catalog Number:	A8300
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BCLAF1 (Phospho-Ser531)
Concentration:	1 mg/ml
Gene Name:	BCLAF1 BTF KIAA0164
Protein Name:	BCLAF1 (Phospho-Ser531)
Other Name:	Bcl-2-associated transcription factor 1 (Btf)
Human Swiss Prot No.:	Q9NYF8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.208% sodium azide.
Specificity:	This antibody detects endogenous levels of BCLAF1 (Phospho-Ser531) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year