

Antibodies Product Datasheet

Product Name:	CAD (Phospho-Thr456) Antibody
Catalog Number:	A8308
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CAD around the phosphorylation site of Thr456. AA range:422-471
Concentration:	1 mg/ml
Gene Name:	CAD
Protein Name:	CAD protein
Other Name:	CAD; CAD protein
Human Swiss Prot No.:	P27708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-CAD (T456) Polyclonal Antibody detects endogenous levels of CAD protein only when phosphorylated at T456.
Storage and Stability:	-20°C/1 year