

Antibodies Product Datasheet

Product Name:	CdGAP (Phospho-Thr789) Antibody
Catalog Number:	A8312
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human CdGAP (phospho Thr789)
Concentration:	1 mg/ml
Gene Name:	ARHGAP31
Protein Name:	Rho GTPase-activating protein 31
Other Name:	ARHGAP31; CDGAP; KIAA1204; Rho GTPase-activating protein 31; Cdc42 GTPase-activating protein
Human Swiss Prot No.:	Q9ULL6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-CdGAP (T789) Polyclonal Antibody detects endogenous levels of CdGAP protein only when phosphorylated at T789.
Storage and Stability:	-20°C/1 year