

Antibodies Product Datasheet

Product Name:	CEP55 (Phospho-Ser425) Antibody
Catalog Number:	A8314
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human CEP55 (phospho Ser425)
Concentration:	1 mg/ml
Gene Name:	CEP55
Protein Name:	Centrosomal protein of 55 kDa
Other Name:	CEP55; C10orf3; URCC6; Centrosomal protein of 55 kDa; Cep55; Up-regulated in colon cancer 6
Human Swiss Prot No.:	Q53EZ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-CEP55 (S425) Polyclonal Antibody detects endogenous levels of CEP55 protein only when phosphorylated at S425.
Storage and Stability:	-20°C/1 year