

Antibodies Product Datasheet

Product Name:	Claspin (Phospho-Ser30) Antibody
Catalog Number:	A8315
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Claspin (Phospho-Ser30)
Concentration:	1 mg/ml
Gene Name:	CLSPN
Protein Name:	Claspin (Phospho-Ser30)
Other Name:	Claspin (hClaspin)
Human Swiss Prot No.:	Q9HAW4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.195% sodium azide.
Specificity:	This antibody detects endogenous levels of Claspin (Phospho-Ser30) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year