

Antibodies Product Datasheet

Product Name:	CRMP-1 (Phospho-Tyr504) Antibody
Catalog Number:	A8323
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CRMP-1 (Phospho-Tyr504)
Concentration:	1 mg/ml
Gene Name:	CRMP1 DPYSL1 ULIP3
Protein Name:	CRMP-1 (Phospho-Tyr504)
Other Name:	Dihydropyrimidinase-related protein 1 (DRP-1) (Collapsin response mediator protein 1) (CRMP-1) (Unc-33-like phosphoprotein 3) (ULIP-3)
Human Swiss Prot No.:	Q14194
Formulation:	Liquid in PBS containing 50% glycerol, and 0.183% sodium azide.
Specificity:	This antibody detects endogenous levels of CRMP-1 (Phospho-Tyr504) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year