

Antibodies Product Datasheet

Product Name:	CUTL1 (Phospho-Ser1215) Antibody
Catalog Number:	A8329
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CUTL1 (Phospho-Ser1215)
Concentration:	1 mg/ml
Gene Name:	CUX1 CUTL1
Protein Name:	CUTL1 (Phospho-Ser1215)
Other Name:	Homeobox protein cut-like 1 (CCAAT displacement protein) (CDP) (Homeobox protein cux-1)
Human Swiss Prot No.:	P39880
Formulation:	Liquid in PBS containing 50% glycerol, and 0.196% sodium azide.
Specificity:	This antibody detects endogenous levels of CUTL1 (Phospho-Ser1215) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year