

Antibodies Product Datasheet

Product Name:	Eps15 (Phospho-Tyr849) Antibody
Catalog Number:	A8342
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Eps15 (phospho Tyr849)
Concentration:	1 mg/ml
Gene Name:	EPS15
Protein Name:	Epidermal growth factor receptor substrate 15
Other Name:	EPS15; AF1P; Epidermal growth factor receptor substrate 15; Protein Eps15; Protein AF-1p
Human Swiss Prot No.:	P42566
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Eps15 (Y849) Polyclonal Antibody detects endogenous levels of Eps15 protein only when phosphorylated at Y849.
Storage and Stability:	-20°C/1 year