

Antibodies Product Datasheet

Product Name:	FKHR (Phospho-Ser249) Antibody
Catalog Number:	A8345
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FKHR (Phospho-Ser249)
Concentration:	1 mg/ml
Gene Name:	FOXO1 FKHR FOXO1A
Protein Name:	FKHR (Phospho-Ser249)
Other Name:	Forkhead box protein O1 (Forkhead box protein O1A) (Forkhead in rhabdomyosarcoma)
Human Swiss Prot No.:	Q12778
Formulation:	Liquid in PBS containing 50% glycerol, and 0.98% sodium azide.
Specificity:	This antibody detects endogenous levels of FKHR (Phospho-Ser249) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year