

Antibodies Product Datasheet

Product Name:	JDP2 (Phospho-Thr148) Antibody
Catalog Number:	A8356
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JDP-2 around the phosphorylation site of Thr148. AA range:114-163
Concentration:	1 mg/ml
Gene Name:	JDP2
Protein Name:	Jun dimerization protein 2
Other Name:	JDP2; Jun dimerization protein 2
Human Swiss Prot No.:	Q8WYK2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-JDP2 (T148) Polyclonal Antibody detects endogenous levels of JDP2 protein only when phosphorylated at T148.
Storage and Stability:	-20°C/1 year