

Antibodies Product Datasheet

Product Name:	Kizuna (Phospho-Thr379) Antibody
Catalog Number:	A8357
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Kizuna (Phospho-Thr379)
Concentration:	1 mg/ml
Gene Name:	PLK1S1 C20orf19 KIZ NCRNA00153 HT013
Protein Name:	Kizuna (Phospho-Thr379)
Other Name:	Centrosomal protein kizuna (Polo-like kinase 1 substrate 1)
Human Swiss Prot No.:	Q2M2Z5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.211% sodium azide.
Specificity:	This antibody detects endogenous levels of Kizuna (Phospho-Thr379) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year