

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                                 |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | LAB (Phospho-Tyr136) Antibody                                                                                                                                                                                                                                   |
| <b>Catalog Number:</b>        | A8359                                                                                                                                                                                                                                                           |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                                                                                                               |
| <b>Applications:</b>          | WB                                                                                                                                                                                                                                                              |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                                                                                                                                                                   |
| <b>Modification:</b>          | Phospho Specific                                                                                                                                                                                                                                                |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                                               |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                                       |
| <b>Immunogen:</b>             | Synthesized peptide derived from human LAB (Phospho-Tyr136)                                                                                                                                                                                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                                         |
| <b>Gene Name:</b>             | LAT2 LAB NTAL WBS15 WBSCR15 WBSCR5 HSPC046                                                                                                                                                                                                                      |
| <b>Protein Name:</b>          | LAB (Phospho-Tyr136)                                                                                                                                                                                                                                            |
| <b>Other Name:</b>            | Linker for activation of T-cells family member 2 (Linker for activation of B-cells) (Membrane-associated adapter molecule) (Non-T-cell activation linker) (Williams-Beuren syndrome chromosomal region 15 protein) (Williams-Beuren syndrome chromosomal region |
| <b>Human Swiss Prot No.:</b>  | Q9GZY6                                                                                                                                                                                                                                                          |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.146% sodium azide.                                                                                                                                                                                                 |
| <b>Specificity:</b>           | This antibody detects endogenous levels of LAB (Phospho-Tyr136) at Human, Mouse, Rat                                                                                                                                                                            |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                                    |