

Antibodies Product Datasheet

Product Name:	MKP-7 (Phospho-Ser446) Antibody
Catalog Number:	A8374
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MKP-7 (Phospho-Ser446)
Concentration:	1 mg/ml
Gene Name:	DUSP16 KIAA1700 MKP7
Protein Name:	MKP-7 (Phospho-Ser446)
Other Name:	Dual specificity protein phosphatase 16 (EC 3.1.3.16) (EC 3.1.3.48) (Mitogen-activated protein kinase phosphatase 7) (MAP kinase phosphatase 7) (MKP-7)
Human Swiss Prot No.:	Q9BY84
Formulation:	Liquid in PBS containing 50% glycerol, and 0.212% sodium azide.
Specificity:	This antibody detects endogenous levels of MKP-7 (Phospho-Ser446) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year