

Antibodies Product Datasheet

Product Name:	MLK3 (Phospho-Ser674) Antibody
Catalog Number:	A8375
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human MLK3 (phospho Ser674)
Concentration:	1 mg/ml
Gene Name:	MAP3K11
Protein Name:	Mitogen-activated protein kinase kinase kinase 11
Other Name:	MAP3K11; MLK3; PTK1; SPRK; Mitogen-activated protein kinase kinase kinase 11; Mixed lineage kinase 3; Src-homology 3 domain-containing proline-rich kinase
Human Swiss Prot No.:	Q16584
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-MLK3 (S674) Polyclonal Antibody detects endogenous levels of MLK3 protein only when phosphorylated at S674.
Storage and Stability:	-20°C/1 year