

Antibodies Product Datasheet

Product Name:	NACA (Phospho-Thr159) Antibody
Catalog Number:	A8380
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NACA (Phospho-Thr159)
Concentration:	1 mg/ml
Gene Name:	NACA HSD48
Protein Name:	NACA (Phospho-Thr159)
Other Name:	Nascent polypeptide-associated complex subunit alpha (NAC-alpha) (Alpha-NAC) (allergen Hom s 2)
Human Swiss Prot No.:	Q13765
Formulation:	Liquid in PBS containing 50% glycerol, and 0.213% sodium azide.
Specificity:	This antibody detects endogenous levels of NACA (Phospho-Thr159) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year