

Antibodies Product Datasheet

Product Name:	NFAT2 (Phospho-Ser172) Antibody
Catalog Number:	A8384
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NFAT2 (Phospho-Ser172)
Concentration:	1 mg/ml
Gene Name:	NFATC1 NFAT2 NFATC
Protein Name:	NFAT2 (Phospho-Ser172)
Other Name:	Nuclear factor of activated T-cells, cytoplasmic 1 (NF-ATc1) (NFATc1) (NFAT transcription complex cytosolic component) (NF-ATc) (NFATc)
Human Swiss Prot No.:	O95644
Formulation:	Liquid in PBS containing 50% glycerol, and 0.92% sodium azide.
Specificity:	This antibody detects endogenous levels of NFAT2 (Phospho-Ser172) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year