

Antibodies Product Datasheet

Product Name:	NFAT5 (Phospho-Ser155) Antibody
Catalog Number:	A8386
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFAT5/TonEBP around the phosphorylation site of Ser155. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	NFAT5
Protein Name:	Nuclear factor of activated T-cells 5
Other Name:	NFAT5; KIAA0827; TONEBP; Nuclear factor of activated T-cells 5; NFAT5; T-cell transcription factor NFAT5; Tonicity-responsive enhancer-binding protein; TonE-binding protein; TonEBP
Human Swiss Prot No.:	O94916
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-NFAT5 (S155) Polyclonal Antibody detects endogenous levels of NFAT5 protein only when phosphorylated at S155.
Storage and Stability:	-20°C/1 year