

Antibodies Product Datasheet

Product Name:	Ack1 (Phospho-Tyr857+Tyr858) Antibody
Catalog Number:	A8455
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Ack1 (Phospho-Tyr857+Tyr858)
Concentration:	1 mg/ml
Gene Name:	TNK2 ACK1
Protein Name:	Ack1 (Phospho-Tyr857+Tyr858)
Other Name:	Activated CDC42 kinase 1 (ACK-1) (EC 2.7.10.2) (EC 2.7.11.1) (Tyrosine kinase non-receptor protein 2)
Human Swiss Prot No.:	Q07912
Formulation:	Liquid in PBS containing 50% glycerol, and 0.190% sodium azide.
Specificity:	This antibody detects endogenous levels of Ack1 (Phospho-Tyr857+Tyr858) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year