

Antibodies Product Datasheet

Product Name:	DAPK1 (Phospho-Ser308) Antibody
Catalog Number:	A9017
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DAPK1 (Phospho-Ser308)
Concentration:	1 mg/ml
Gene Name:	DAPK1 DAPK
Protein Name:	DAPK1 (Phospho-Ser308)
Other Name:	Death-associated protein kinase 1 (DAP kinase 1) (EC 2.7.11.1)
Human Swiss Prot No.:	P53355
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of DAPK1 (Phospho-Ser308) at Human:S308, Mouse:S308
Storage and Stability:	-20°C/1 year