

Antibodies Product Datasheet

Product Name:	EPHA1 (Phospho-Tyr605) Antibody
Catalog Number:	A9020
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHA1 (Phospho-Tyr605)
Concentration:	1 mg/ml
Gene Name:	EPHA1 EPH EPHT EPHT1
Protein Name:	EPHA1 (Phospho-Tyr605)
Other Name:	Ephrin type-A receptor 1 (hEpha1) (EC 2.7.10.1) (EPH tyrosine kinase) (EPH tyrosine kinase 1) (Erythropoietin-producing hepatoma receptor) (Tyrosine-protein kinase receptor EPH)
Human Swiss Prot No.:	P21709
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of EPHA1 (Phospho-Tyr605) at Human:Y605, Mouse:Y606
Storage and Stability:	-20°C/1 year