

Antibodies Product Datasheet

Product Name:	EPHB4 (Phospho-Tyr987) Antibody
Catalog Number:	A9033
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHB4 (Phospho-Tyr987)
Concentration:	1 mg/ml
Gene Name:	EPHB4 HTK MYK1 TYRO11
Protein Name:	EPHB4 (Phospho-Tyr987)
Other Name:	Ephrin type-B receptor 4 (EC 2.7.10.1) (Hepatoma transmembrane kinase) (Tyrosine-protein kinase TYRO11)
Human Swiss Prot No.:	P54760
Formulation:	Liquid in PBS containing 50% glycerol, and 0.109% sodium azide.
Specificity:	This antibody detects endogenous levels of EPHB4 (Phospho-Tyr987) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year