

Antibodies Product Datasheet

Product Name:	Mre11 (Phospho-Ser676) Antibody
Catalog Number:	A9047
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Mre11 (Ser676)
Concentration:	1 mg/ml
Gene Name:	MRE11A HNGS1 MRE11
Protein Name:	Mre11 (Ser676)
Other Name:	Double-strand break repair protein MRE11A (Meiotic recombination 11 homolog 1) (MRE11 homolog 1) (Meiotic recombination 11 homolog A) (MRE11 homolog A)
Human Swiss Prot No.:	P49959
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mre11 (phospho-Ser676)
Storage and Stability:	-20°C/1 year