

Antibodies Product Datasheet

Product Name:	Bamacan (Phospho-Ser1083) Antibody
Catalog Number:	A9069
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Bamacan (Phospho-Ser1083)
Concentration:	1 mg/ml
Gene Name:	SMC3 BAM BMH CSPG6 SMC3L1
Protein Name:	Bamacan (Phospho-Ser1083)
Other Name:	Structural maintenance of chromosomes protein 3 (SMC protein 3) (SMC-3) (Basement membrane-associated chondroitin proteoglycan) (Bamacan) (Chondroitin sulfate proteoglycan 6) (Chromosome-associated polypeptide) (hCAP)
Human Swiss Prot No.:	Q9UQE7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of Bamacan (Phospho-Ser1083) at Human:S1083, Mouse:S1083, Rat:S1083
Storage and Stability:	-20°C/1 year