

Antibodies Product Datasheet

Product Name:	EIF5 (Phospho-Ser389+Ser390) Antibody
Catalog Number:	A9081
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF5 (Phospho-Ser389+Ser390)
Concentration:	1 mg/ml
Gene Name:	EIF5
Protein Name:	EIF5 (Phospho-Ser389+Ser390)
Other Name:	Eukaryotic translation initiation factor 5 (eIF-5)
Human Swiss Prot No.:	P55010
Formulation:	Liquid in PBS containing 50% glycerol, and 0.200% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF5 (Phospho-Ser389+Ser390) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year