

Antibodies Product Datasheet

Product Name:	EIF6 (Phospho-Ser235) Antibody
Catalog Number:	A9082
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF6 (Phospho-Ser235)
Concentration:	1 mg/ml
Gene Name:	EIF6 EIF3A ITGB4BP OK/SW-cl.27
Protein Name:	EIF6 (Phospho-Ser235)
Other Name:	Eukaryotic translation initiation factor 6 (eIF-6) (B(2)GCN homolog) (B4 integrin interactor) (CAB) (p27(BBP))
Human Swiss Prot No.:	P56537
Formulation:	Liquid in PBS containing 50% glycerol, and 0.155% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF6 (Phospho-Ser235) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year